+tucky

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OPERATING INSTRUCTIONS e-health thermometer model TUCKY-21

INTRODUCTION

The Tucky thermometer has been designed to continuously and remotely measure the child's temperature and position.

It is recommended for children from birth to 5 years of age.

It is positioned under the armpit with a dedicated, double-sided adhesive, and sends the information gathered to a smartphone application.

The temperature reading displayed is colour coded, indicating whether the temperature is normal or not.

Correspondence table for colour/axillary temperature:

Green Normal temperature	Orange High temperature	Red Very high temperature	Grey Temperature not representative
34.7°C to 37.3°C	37.3°C to 38.4°C (default value)	higher than 38.4°C (default value)	E.g. arm open*

NB: Rectal temperature is usually 0.6°C higher than axillary temperature.

The thermometer is adjusted in direct mode, i.e. it takes the actual axillary temperature, without making adjustment calculations. You should therefore wait about fifteen minutes after fitting Tucky in order for the temperature to stabilise.

* Axillary temperature varies according to the individual person (age, sex, type, skin colour), the time of day and external factors (ambient temperature, clothing etc.) [1]. The measurements are least disturbed by environmental conditions when the child's arm is closed. Tucky ensures continuous monitoring. **It is normal to observe variations in temperature caused by changes in environmental conditions or the position of the arm.** However, the trend indicated by the curve remains significant. The measurement is considered not to be representative when it is displayed in grey.

The default "very high temperature" value is set at 38.4°C and this value can be adjusted by the user based on their knowledge of their child's physiological data.

In case of Bluetooth disconnection, the data is no longer displayed. The connection is re-established automatically when the communication is possible. As with all thermometers, it is essential to use the right technique to obtain accurate temperatures. Therefore, please read the instructions carefully and in full before using Tucky.

WARNINGS AND PRECAUTIONS FOR USE

This product is solely for domestic use and child-minding.

The Tucky medical device is not designed for diagnosis, it's a monitoring tool. The use of this medical device must under no circumstances replace supervision of one's child or consulting a doctor. If unusual signs or symptoms appear in the child (vomiting, diarrhoea, dehydration etc.), even in the absence of a temperature, parents must consult a doctor.

Consult your doctor if the temperature displayed is high. A persistent high fever requires medical treatment, particularly in young children.

The device must not be used beyond the intended life cycle. Similarly, it is not recommended for the child to wear Tucky for extended periods: it is recommended that it is not worn for longer than 24 hours.

Tucky may be used at an ambient temperature of between +15°C and +35°C and stored at temperatures of

between -15°C and +55°C: keep this thermometer away from direct sunlight, fire and any contamination. Do not use this thermometer to interpret hypothermia temperatures.

It must not be submerged in water: take Tucky off the child before a bath.

It must be used with the strip-cap **1** closed.

When it is charging, the device must under no circumstances be placed on your child. Similarly, this device is intended to be used by adults and must under no circumstances be left within reach of children when it is not being used.

To charge the Tucky, you must connect a cable (see "Technical characteristics") to a separate 5V power supply (USB port). You must use the cable provided or an equivalent bearing CE marking.

Tucky and its adhesives are designed with biocompatible medical materials. However, do not use this product if your child has already been diagnosed with an allergy to elastomer or silicon type plastics.

The product has been designed to fit snugly to the curves of the body. Do not bend it other than for its intended use.

Before using Tucky, make sure the time is properly set on your phone or tablet. The information issued by the app is based on the time entered on your telephone or tablet.

It is advised to apply Tucky to clean skin, and if applicable, to wait for the child's skin to have absorbed cream or oil after application. Do not attach to irritated or wounded skin.

Do not modify this device without the manufacturer's authorisation.

Tucky is provided with the following accessories: a box of adhesives (ref. ADH-21) and a USB lead, which must not, under any circumstances, be left within reach of children.

The device must be applied to the child's skin with the adhesives provided in the packaging. Under no circumstances use other products to do this. The adhesives must be stored at a temperature of between +10°C and +27°C, away from sunlight and in a dry place in order to preserve their properties.

Adhesives are available on the www.e-takescare.com website and in most shops selling Tucky. They cannot be repositioned and are single-use only.

The manufacturing date and SN serial number can be accessed via a barcode (QRcode or Datamatrix).

Apps for reading barcodes are available on the Apple App Store and the Google Play Store.

This information is on the label, situated under the strip-cap **1**.

The manufacturer reserves the right to modify the product's technical characteristics without prior notice.

TRANSMISSION OF INFORMATION

2 scenarios:

Direct reading: the smartphone or tablet is connected to Tucky via Bluetooth (to read data from around 10 m away).

Distance reading: to read data from further away than 10 m, a second device (smartphone or tablet with the app installed) must be placed within Bluetooth reach of Tucky in order to relay the data to the user via the internet. The data is stored on an accredited healthcare host.

HOW TO USE TUCKY

First use

Charge Tucky for at least 1 h 30 mins, using the charger lead provided. Plug the USB micro connector into Tucky and the USB connector to a computer, a tablet or a USB power adapter. 10 minutes of charging provides 8 hours of use approximatively.

Download the free «Tucky» app from the Apple App Store or Google Play Store (check the "Technical characteristics" paragraph beforehand to see if it is compatible with your phone).

Activate Bluetooth.

Open the app: create a user account following the instructions on the screen.

Then, at each use

- Clean Tucky (as explained in the «maintenance and cleaning» paragraph).
- Unstick a double-sided adhesive from its backing and apply the sticky side to the side with the metal sensor, taking care to leave the sensor uncovered.
- Remove the protective film from the adhesive placed on the Tucky. Place Tucky under the armpit (see "How to position Tucky under the child's armpit" paragraph).
- Press the on/off button (the LED will blink).
- If Tucky is already linked to a profile, the connection is made automatically. If not, select a profile for your child. Select the Menu icon. Choose Tucky in the list of devices detected.


Tucky is now operational.

To turn off Tucky, press the on/off button (the LED will light up for 2 seconds). Otherwise, Tucky will switch off automatically 30 minutes after being removed from the child.

How to position Tucky under the child's armpit

In order to guarantee that the temperature is taken as accurately as possible, Tucky must be placed according to these 3 points:

- under the right armpit, as close as possible to the armpit,
- with the temperature sensor  positioned on the child's skin,
- with the micro-USB port positioned downwards.

Maintenance and cleaning

Use a cloth or a bit of cotton wool soaked in isopropyl or ethyl alcohol to clean Tucky.

Do not use abrasive cleaning products.

This thermometer must not be washed in the dishwasher or with boiling water. If it is washed in either of these ways, Tucky will no longer work and this will void the warranty.

Alerts and notifications

The application signals an alert in the following cases:

- temperature > 38.4°C for longer than 5 minutes,
- low battery,
- Tucky detached,
- loss of long connection,
- malfunction (error).

Technical characteristics

Product type	Medical thermometer in direct mode
Reference	Tucky-21
Medical device	Classe IIa
Normal condition of use	Ambient operating temperature: +15°C to 35°C Technical accuracy cannot be guaranteed if the device is used outside the temperature and humidity ranges indicated
Storage/transport conditions	Tucky-21W: Temperature: -15°C to 55°C ± 2°C. Humidity: 15% to 95% ± 5% RH ADH-21: Temperature: +10°C to 27°C ± 2°C. Humidity: 40% to 60% ± 5% RH
Atmospheric pressure	Can operate at an atmospheric pressure which does not exceed 1 (760 à 1060 hPa)
Rate of humidity	Relative humidity range of 15% to 85% (without condensation)
Battery	Device with internal power supply: lithium battery. Capacity: 15mAh
Range of temperature displayed	+24°C to +44°C
Display resolution	0.1°C
Accuracy	Accuracy for the range of temperatures displayed (maximum error in the laboratory): ±0.2°C from 32°C to 42 °C then +/- 0.3°C outside this range
Technologie Bluetooth	4.0
Compatibility	Compatibility: from iOS 7.3 and Android 4.3
USB Cable	USB 2.0 Cable Type A Male to Type B Male
Charge voltage	5V
Protection index	IP22: Protected against solid objects larger than 12 mm. Protected against vertically falling drops of water up to 15°C
Lifetime of Tucky	500 charging cycles
Lifetime of ADH-21 adhesives	2 years
Software classification	Class B

ESSAIS ELECTROMAGNETIC EMISSION TESTING	COMPLIANCE	ENVIRONMENT
RF Emissions (CISPR 11)	Group 1	The RF emissions are very low and are not likely to cause interference in nearby electronic equipment.
RF Emissions (CISPR 11)	Class B	Compliant
Harmonic current emissions	Not applicable	Battery power
Voltage fluctuation, flicker emissions	Not applicable	Battery power

IMMUNITY TESTING	TESTING LEVEL ACCORDING TO CEI 60601	LEVEL OF COMPLIANCE
Electrostatic discharge (ESD) CEI 61000-4-2	± 8 kV at contact ± 2,4,8 and 15 kV in the air	Compliant
RF radiated IEC 60601-4-3	10V/m de 80 MHz à 2,5 GHz	Compliant
Magnetic field at the network frequency IEC 60601-4-8	30 A/m	Compliant

Warranty : please refer to the general sales terms, available on the website www.e-takescare.com

Not intended to be sterilised. Subject to modifications without notice.

Inspection, repair and calibration

At each battery charge, an auto-diagnosis is performed to check Tucky's operation. When Tucky no longer works, an error message is sent to the app.

This device has been designed and manufactured to last. It is nevertheless recommended to have it inspected once a year by e-TakesCare, to guarantee its operation and accuracy.

Checking accuracy is not free of charge. If applicable, we encourage you to make contact with e-takescare.

Memory and connection

Before using Tucky, you must ensure that there is enough memory on the smartphone or tablet to store the data if the Internet disconnects. As Tucky's memory is limited, it is recommended to avoid Bluetooth disconnection phases as far as possible. If you encounter problems with the Bluetooth connection: check that Tucky is within reach of the device you wish to sync with, and check that Bluetooth is enabled on your device.

Regulations This device is compliant with the following standards:

NF EN ISO 80601-2-56: 2012 Electromedical device - Part 2-56: Specific requirements relating to basic safety and essential performance of medical thermometers for measuring body temperature.

NF EN ISO 14971: 2013 Medical devices – Application of the management of risks to medical devices.

NF EN 62366: 2008 + A1: 2015 Medical devices – Application of engineering of aptitude for use with medical devices.

NF EN 60601-1: 2006 + A1 Electromedical devices - Part 1: General requirements for basic safety and essential performance + Appendix I.

NF EN 60601-1-2: 2007 Electromedical devices - Part 1-2: General requirements for basic safety and essential performance + Collateral standard: Electromagnetic compatibility.

NF EN 60601-1-1: 2010 Electromedical devices - Part 1-1: General requirements for basic safety and essential performance + Collateral standard: Requirements for electromedical devices and electromedical systems used in the home care environment.

NF EN 62304: 2006 Medical device software. Software life cycle process.

NF EN ISO 15223-1: 2012 Symbols to use with labels, labelling and the information to provide relating to medical devices - Part 1: General requirements.

NF EN 1041: 2008 + A1: 2013 Information provided by the manufacturer of medical devices.

ETSI EN 300 328 V1-9-1: 2012-06 Bluetooth

IEC 62133: Battery certificate

EN 71-1: Toy safety: «Testing size» part

Portable of mobile RF communication devices may interfere with electro-medical devices.

[1] KELLY, Greg. Body temperature variability (Part 1): A review of the history of body temperature and its variability due to site selection, biological rhythms, fitness, and aging. *Alternative medicine review*, 2006, vol. 11, no 4, p. 278.

After-sales service For any technical problems, please visit the www.e-takescare.com website

MEANING OF SYMBOLS

CE labelling given in 2016

Non-ionising radiation

Direct current

Follow the instructions for use

Device containing BF type applied parts

Manufacturer: e-TakesCare
35, rue des Chantiers
F - 78000 Versailles

Storage temperature
-15°C to 55°C

The device, its accessories and its packaging must be recycled appropriately at the end of their use. Please observe local regulations and rules

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