

+LUCKY



SMART WEARABLE THERMOMETER


Safe fever monitoring

Continuous monitoring

In the form of a soft, comfortable and supple **patch**, Tucky is placed under the armpit and allows to monitor **evolution of body temperature**.

Real time temperature curve is displayed on a remote smartphone app.

An **alert** is raised when fever exceeds a threshold temperature (audible and visual notification).



Without distance limit

Fever can be monitored either at home (5-10m via Bluetooth) or remotely via internet on a 2nd smartphone (secured cloud data synchronization).



A simple and intuitive app

Multi profiles, multi devices, curve / timeline views , data history,, health book...

A secured medical device



Class IIA Medical Device



Patented Technology

Ultra low wave emission



Secured data

Certified Healthcare Cloud

The most complete thermometer of the market

Thermometer	tucky	Ear / temporal	Traditional
			
24h continuous monitoring	✓	x	x
Remote tracking	✓	x	x
Fever alert	✓	✓	x
Temperature color code	✓	✓	x
Data history	✓	✓	x
Position monitoring	✓	x	x

- Flexible, comfortable, hypoallergenic
- Medical certification, class Ila, accuracy 0.1°C
- Ultra low wave emission (x1000 lower than a smartphone)
- Rechargeable (autonomy > 5 day)
- For babies (>29 days), children and adults.



+ product

Tucky is also a **position monitor**. It warns when newborn is no longer in the recommended sleeping position.

Tucky box & adhesive refill kit

Tucky box includes:

- 1 wearable thermometer
- 5 double-faces adhesives
- 1 USB recharge cable
- 1 manual

Adhesive refill kit: x15 adhesives (ref: ADH-21-P1)



Model TUCKY-21W – Class Ila medical device
 Manufactured by e-TakesCare SAS, 35 rue des Chantiers, 78000 Versailles, France
 More information: www.e-takescare.com

